

Tristel Wipes System

The Tristel Wipes System is a practical and highly effective way to decontaminate and disinfect non-lumened medical devices, including those that cannot be fully immersed in liquid disinfectant or sterilants and cannot be sterilised by heat. Tristel Wipes System is not intended for use on critical medical devices which must be sterilised prior to use on a patient.

A three part system

- Pre-Clean wipe
- Sporidical wipe
- Rinse wipe

Features

- Sporidical wipes kill all organisms on a pre-cleaned surface in 30 seconds
- Sporidical wipe is a Class IIb Medical Device
- All wipes are individually packed in sachets
 - Pre-clean Wipes
 - Sporidical Wipes
 - Rinse Wipes
- Each Sporidical wipe has a unique traceability label enabling Quality Audit Trail
- Rinse wipe is sterile packed
- Easy to use – time saving
- More economical than ‘sheaths’

Applications

- Nasendoscopes
- Transoesophageal echo (TOE) probes
- Trans-vaginal and trans-rectal ultrasound probes
- G.I. High resolution manometry catheters
- Ophthalmic devices

Patent no: GB 2404 337B & WO 2005/011756A1



Code: TRIO25

1 carton: 25 of each wipe + 1 Activator Foam

Code: TRIO50

1 carton: 50 of each wipe + 1 Activator Foam

Code: PRE/W - TSW₁ - RIN/W

1 carton of each type: 50 wipe sachets


Tristel Wipes System

The decontamination process involves three stages:

1. Cleaning

2. Disinfection

3. Rinsing



Tristel Pre-Clean Wipe

The first step in the decontamination process of medical devices is the thorough cleaning of the surface to remove soil and organic matter prior to high-level disinfection/sterilisation.

The Tristel Pre-Clean Wipe is impregnated with a low-foaming surfactant system combined with triple enzymes, producing ultra-low surface tension for rapid cleaning of any hard surface.

The Tristel Pre-Clean Wipe is a Class I medical device carrying CE mark in accordance with the European Medical Device Directive 93/42/EEC.

Tristel Sporicidal Wipe

CE marked as a Class IIb Medical Device, Tristel Sporicidal Wipes are the central part of the Tristel Wipes System and have been designed specifically for high-level sporicidal disinfection ENT scopes, ultrasound probes and other non-lumened medical devices. The wipe, once activated with the foam, is sporicidal, mycobactericidal, bactericidal, virucidal and fungicidal within a 30 second contact time.

Tristel Rinse Wipe

The final step in the decontamination process is the rinsing of the surface that has been treated by a chemical biocide. The Tristel Rinse Wipe is a sterile-packed and is impregnated with de-ionised water. It is used to remove chemical residues from a surface that has been decontaminated with a Tristel Sporicidal Wipe.

The Tristel Rinse Wipe is a Class I Sterile medical device carrying CE mark in accordance with the European Medical Device Directive 93/42/EEC.

Tristel Wipe Quality Audit Trail

Enables the user to meet the directive of the Chief Medical Officer in “Winning Ways” (Working together to reduce Hospital Acquired Infection), December 2003, in which it states that staff will ensure that there is an audit trail for each recycled item.