

# Tristel

## Fuse for Labs

Sporicidal disinfectant



### SPORICIDAL DISINFECTANT FOR USE IN LABS:

- **In-vitro diagnostic medical devices**
- **Discard pots**
- **Centrifuges**
- **Spillages**
- **Isolators**
- **Work surfaces**
- **Benches**
- **Floors and walls**

#### Each sachet contains:

- 50ml Base solution (solution of citric acid, with preservatives, surfactant and corrosion inhibitors in demineralised water)
- 50ml Activator solution (sodium chlorite solution in demineralised water)

This makes 100ml chlorine dioxide in aqueous solution for dilution with water.

### Chlorine dioxide is effective against:

Spores	<i>Bacillus cereus</i> <i>Bacillus subtilis</i> <i>Clostridium difficile</i>
Mycobacteria	<i>Avium</i> <i>Terrae</i> <i>Tuberculosis</i>
Viruses	<i>Adenovirus</i> <i>Hepatitis</i> <i>Herpes</i> <i>HIV</i> <i>Orthopox</i> <i>Polio</i>
Fungi	<i>Aspergillus niger</i> <i>Candida albicans</i>
Bacteria	<i>Enterococcus faecium</i> <i>VRE</i> <i>Enterococcus hirae</i> <i>Escherichia coli</i> <i>Pseudomonas aeruginosa</i> <i>GRP</i> <i>Staphylococcus aureus</i> <i>MRSA</i>

# Tristel

## How to use Tristel Fuse for Labs



1. Take one sachet to produce two litres of  $\text{ClO}_2$  solution



2. Squeeze one side of the sachet to burst contents through centre seal. Contents will start to turn yellow




3. Move contents from side to side for 30 seconds



4. Cut corner of sachet



5. Pour contents into two litres of water



6. Sporidical in five minutes



### Manufactured by:

Tristel Solutions Limited  
Cams, CB8 7NY, UK  
Telephone: +44 (0)1638 721500  
[www.tristel.com](http://www.tristel.com)

