



NEW RELEASE
22 December 2005

Tristel Duo sporicidal cleaner effective against *C. difficile*

Current concern about the risk of *Clostridium difficile* infection in the hospital environment highlights the importance of having effective cleaning and sterilisation products available for immediate use. Tristel's new Duo wipes system offers convenient, highly effective bactericidal, fungicidal and sporicidal cleaning for use on all general hospital hard surfaces and non-invasive medical devices. It is effective against a wide range of organisms including *Clostridium difficile*.

Effective, convenient and easy to use, Tristel Duo is designed to provide a readily available means of rapid decontamination and is active against organisms including MRSA, *Clostridium difficile*, *Bacillus subtilis*, *Mycobacterium tuberculosis*, *Acinetobacter baumannii* and *Norovirus*. It has applications in diverse areas, including for work surfaces, dressing trolleys, beds and mattress covers, instrument tables, monitor cases, control panels and keyboards. The time required to kill all the organisms that Tristel Duo comes into contact with is just 30 seconds.

Tristel Duo is a two part system. It comprises the Duo bactericidal and virucidal wipe, which can be boosted by the application of the Tristel Duo sporicidal foam. This delivers chlorine dioxide directly to the surface, as and when required

For full effectiveness, chlorine dioxide must be generated immediately before use from two chemical components. Tristel has achieved this using the Duo foamer, a novel two-compartment dispenser. Base and activator solutions are held in separate compartments of the Duo foamer and are mixed only when the pump is depressed. Two jets of foam are dispensed to the surface where they are mixed and dispersed with the Duo wipe.

For more details of the complete range of Tristel products visit www.tristel.com , call +44 1638 721500 or email : mail@tristel.com

Ends

Illustration, notes and contact details follow...

Tristel Duo sporicidal cleaner effective against *C. difficile*.../2

High resolution image available from Kapler Communications
trish@kapleronline.com Please quote reference TRI/PRL/065



About Tristel Solutions Limited

The cornerstone of the Tristel product range is a patented chlorine dioxide chemistry used to sterilise medical, dental and veterinary instruments. In 1998, the business that became Tristel acquired all worldwide rights to this chemistry. The principal operating subsidiary of the group was established at this time to exploit the UK market and Tristel has become the leading supplier of liquid chemical sterilising solutions to UK hospitals.

The recently introduced Tristel Wipe is the world's first rapid action sporicidal wipe. It incorporates Tristel's patented chlorine dioxide chemistry and can kill all organisms on a pre-cleaned surface, from which soil and organic matter have been removed, with a contact time of only 30 seconds. The Tristel Generator and Tristel Duo are the latest addition to the company's range of products and services.

Tristel Solutions Limited was admitted to trading on AIM on 1 June 2005. Its stock exchange symbol is TSTL.

www.tristel.com

Press information

Trish Appleton
Kapler Communications
Suite 2, Cressner House
12 Huntingdon Street
St Neots
Cambs UK
PE19 1BD

t: +44 (0) 1480 471117
f: +44 (0) 1480 471118
trish@kapleronline.com

Sales enquiries

Polly Oates
Tristel Solutions Limited
Lynx Business Park
Fordham Road
Snailwell
Cambs UK
CB8 7NY

t: +44 (0) 1638 721500
f: +44 (0) 1638 721911
mail@tristel.com
www.tristel.com