

# Tristel Jet™

## Liquid for Surfaces

### Microbiological Efficacy Summary Australia and New Zealand

Sporicidal	Test norm	Test type	Organism	Laboratory	Clean or dirty	Contact time
	EN 14561	Carrier	<i>Bacillus subtilis</i>	HygieneNord	Clean 1	2 min.
			<i>Clostridium difficile</i>	HygieneNord	Clean 1	

Mycobactericidal	Test norm	Test type	Organism	Laboratory	Clean or dirty	Contact time
	EN 14563	Carrier	<i>Mycobacterium avium</i>	HygieneNord	Clean 1	60 sec.
	EN 14348	Suspension	<i>Mycobacterium terrae</i>	MGS	Clean 1	

Virucidal	Test norm	Test type	Organism	Laboratory	Clean or dirty	Contact time
	EN 14476	Suspension	Poliovirus Type 1	Dr Brill & Dr Steinmann	Dirty 1	30 sec.
			Adenovirus Type 5	Dr Brill & Dr Steinmann	Dirty 1	
			Human Herpesvirus Type 1	Dr Brill & Dr Steinmann	Dirty 1	30 sec.

Yeasticidal	Test norm	Test type	Organism	Laboratory	Clean or dirty	Contact time
	EN 14562	Carrier	<i>Candida albicans</i>	MGS	Clean 1	60 sec.

Bactericidal	Test norm	Test type	Organism	Laboratory	Clean or dirty	Contact time
	EN 14561	Carrier	<i>Pseudomonas aeruginosa</i>	MGS	Clean 1	60 sec.
			<i>Enterococcus hirae</i>	MGS	Clean 1	
			<i>Staphylococcus aureus</i>	MGS	Clean 1	
	EN 13727	Surface	Carbapenemase-resistant Enterobacteriaceae (CRE) <i>Klebsiella pneumoniae</i>	MGS	Clean 1	
			Vancomycin-resistant Enterococci (VRE) <i>Enterococcus faecium</i>	MGS	Clean 1	
			Multi-drug resistant <i>Acinetobacter baumannii</i> (MRDAB)	MGS	Clean 1	
			Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA)	MGS	Clean 1	

#### Clean/Dirty Conditions Key

Clean 1: 0.3 g/L bovine albumin

Dirty 1: 10% Fetal Calf Serum



Created by: Tristel Solutions Limited, Lynx Business Park, Cambs, UK, CB8 7NY  
T +44 (0) 1638 721500 - E mail@tristel.com - W www.tristel.com

Australian sponsor: Tristel Pty Ltd, 40/328 Reserve Road, Cheltenham, VIC 3192  
T 1300 680 898 - F +61 (0)3 9533 6193 - E mail-au@tristel.com

New Zealand sponsor: Tristel New Zealand Limited, P.O Box 15451, Tauranga, New Zealand  
T +64 (0) 7 5773560 - E info@tristel.co.nz

Copyright © Tristel Solutions  
Mkt-Mic-431-1  
August 2018

For Tristel patent information please visit: <http://www.our-patents.info/tristel>