

# DISTEL

— c —



One bottle of DISTEL concentrate solution generates 2000 surface disinfectant wipes.

DISTEL is intended for the cleaning and disinfection of non-invasive medical devices and hard surfaces within the hospital environment.

DISTEL is effective against bacteria, fungi and viruses including coronaviruses. It is also effective against blood borne viruses (HIV/ Hepatitis B) and destroys DNA and RNA.

**ELIMINATE SINGLE-USE PLASTIC WIPES**



**MANUFACTURED IN THE UK**




**REUSABLE BOTTLES**




**USE THE WIPE OF YOUR CHOICE**

## SAVING SPACE, REDUCING WASTE



10 x 200 IMPREGNATED WIPES PACKS




=

1 x 1 LITRE DISTEL CONCENTRATE

DISTEL can be used with any wipe to create freshly made surface disinfectant wipes on the go. Choose a sustainably sourced wipe to reduce single-use plastic waste.

## ENVIRONMENTALLY CONSCIOUS


Designed to be repurposed and reused, DISTEL dosing bottles keep the environment in mind. One DISTEL concentrate bottle delivers 10 x 1 litre bottles of working solution which equates to 200 x 5ml doses per bottle.



**DOSE**



**WIPE**



**REPEAT\***


cache™

TRISTEL'S CLEANING AND DISINFECTION BRAND FOR SURFACES


**DISTEL IS INTENDED FOR THE CLEANING AND DISINFECTION OF NON-INVASIVE MEDICAL DEVICES AND HARD SURFACES WITHIN THE HOSPITAL ENVIRONMENT SUCH AS:**



**WORKSTATIONS  
(MONITORS & KEYBOARDS)**




**TABLE TOPS**




**TROLLEYS**



**HOISTS**




**HANDRAILS**



x5      x1      x6

**STARTER PACK  
CCH050101**


1 x 1 Litre Bottle of Concentrate  
5 x 1 Litre Dosing Bottle  
6 x 5ml Dosing Heads  
1 x Concentrate Bottle Label  
6 x Working Solution Bottle Labels



x6      x6

**CONCENTRATE PACK  
CCH050201**

6 x 1 Litre Bottle of Concentrate  
6 x 5ml Dosing Heads  
6 x Concentrate Bottle Label  
6 x Working Solution Bottle Labels



x6      x6

**DOSING PACK  
CCH050301**

6 x 1 Litre Dosing Bottle  
6 x 5ml Dosing Heads  
6 x Working Solution Bottle Labels

**AVAILABLE TO  
ORDER NOW!**

DISTEL Medical Surface Disinfectant is CE marked as a Class IIa Medical Device.

\*Dose onto a new wipe. For full instructions please see the DISTEL User Guide.